

# REDWOOD VISCOMETER

**Ex.No :**

**Date:**

**Aim :**

To determine the kinematic viscosity and absolute viscosity of the given lubricating oil at different temperatures using Redwood Viscometer

## **Apparatus required :**

Redwood Viscometer  
Thermometer 0-100°C (2 Nos)  
Stop watch  
50 ml standard narrow necked flask  
Given Sample of oil


## **Description :**

The redwood viscometer consist of vertical cylindrical oil cup with an orifice in the centre of its base . The orifice can be closed by a ball . A hook pointing upward serve as a guide mark for filling the oil . The cylindrical cup is surrounded by the water bath . The water bath maintain the temperature of the oil to be tested at constant temperature . The oil is heated by heating the water bath by means of an immersed electric heater in the water bath , The provision is made for stirring the water , to maintain the uniform temperature in the water bath and to place the thermometer to record the temperature of oil and water bath . The cylinder is 47.625mm in diameter and 88.90mm deep . The orifice is 1.70mm in diameter and 12mm in length , This viscometer is used to determine the kinematic viscosity of the oil. From the kinematic viscosity the dynamic viscosity is determined .

## **Theory and Definition :**

Viscosity is the property of fluid . It is defined as “The internal resistance offered by the fluid to the movement of one layer of fluid over an adjacent layer ‘ . It is due to the Cohesion between the molecules of the fluid . The fluid which obey the Newton law of Viscosity are called as Newtonian fluid .

The dynamic viscosity of fluid is defined as the shear required to produce unit rate of angular deformation .



### Observation and tabulation

- (1) Room temperature  $T_R = \dots\dots\dots$  °C  
 (2) Density of oil at room temperature =  $\dots\dots\dots$  gm/cm<sup>2</sup>

S.No	Temperature of oil °C	Time taken to fill 50ml flask in 'Sec'	Kinematic Viscosity in 'Centi Stokes'	Density in gm/cc	Dynamic (or) Absolute viscosity 'Centi Poise'
1					
2					
3					
4					
5					
6					

**Formulae used :**

$$\text{Kinematic Viscosity} = \frac{A}{t} - \frac{B}{t} \quad \begin{array}{l} \text{in stokes or} \\ \text{in centi stokes} \end{array}$$

$$A = 0.0026$$

$$B = 1.72$$

$$A = 0.26$$

$$B = 172$$

t = Saybolt second

Density of oil at particular temperature  $\rho_t$

$$\rho_t = \rho_R - 0.00065 (T - T_R)$$

T = Temperature at which the density is required

T<sub>R</sub> = Room Temperature

$\rho_R$  = Density of oil at room temperature in gm / cm<sup>3</sup>  
= 0.84 (or) 0.85 gm/cm<sup>3</sup>

i.e.,  $\tau = \mu \frac{du}{dy}$

$$\mu = \frac{\tau}{du/dy}$$

where  $\mu$  = Co-efficient of viscosity (or)  
Dynamic viscosity (or)  
Absolute viscosity

$\tau$  = Shear stress

du = Angular deformation (velocity gradient )

The unit of dynamic viscosity in SI system is

$$\frac{\text{N} - \text{Sec}}{\text{m}^2} \quad (\text{or}) \quad \frac{\text{kg}}{\text{m} - \text{sec}} \quad (\text{or}) \quad \text{poise}$$

In metric system :

$$\frac{\text{dynes} - \text{Sec}}{\text{m}^2} \quad \text{or} \quad \frac{\text{gm}}{\text{cm} - \text{Sec}}$$

$$\frac{1 \text{ N-S}}{\text{m}^2} = 10 \text{ Poise}$$

The kinematic viscosity of the fluid is defined as the ratio of the dynamic viscosity to the mass density of the fluid. Its symbol is 'ν'

$$\nu = \frac{\mu}{\rho} ; \quad \rho = \text{mass density of oil}$$


The unit of kinematic viscosity

In SI system :  $\text{m}^2/\text{sec}$

In metric system =  $\frac{\text{cm}^2}{\text{sec}}$  or Stokes

$$1 \frac{\text{m}^2}{\text{sec}} = 10^4 \frac{\text{cm}^2}{\text{sec}} = 10^4 \text{ stokes}$$

One hundred<sup>th</sup> part of stoke is called Centi Stoke .




## Procedure :

- (1) Clean the cylindrical oil cup and ensure the orifice tube is free from dirt .
- (2) Close the orifice with ball valve.
- (3) Place the 50 ml flask below the opening of the Orifice .
- (4) Fill the oil in the cylindrical oil cup upto the mark in the cup .
- (5) Fill the water in the water bath.
- (6) Insert the thermometers in their respective places to measure the oil and water bath temperatures.
- (7) Heat the by heating the water bath, Stirred the water bath and maintain the uniform temperature .
- (8) At particular temperature lift the bal valve and collect the oil in the 50 ml flask and note the time taken in seconds for the collecting 50 ml of oil . A stop watch is used measure the time taken . This time is called Redwood seconds .
- (9) Increase the temperature and repeat the procedure '8' and note down the Redwood seconds for different temperatures .

## Graph :

The following graph has to be drawn

- (1)Temperature Vs Redwood seconds
- (2)Temperature Vs Kinematic Viscosity
- (3)Temperature Vs Dynamic Viscosity



## Result :

The kinematic and dynamic viscosity of given oil at different temperatures were determined .